

SNIJDERS CARE manufactures and supplies internal transport systems.

In addition to the extensive range of standard products covered in this brochure, we can meet custom requirements in size and finish, and full customized solutions for new or existing departments.

Our unchanging priority is to deliver the best possible end product based on the user's experience and knowledge of the market.

*Specifications are subject to change.

CONTENTS

PAGE

STAINLESS STEEL LAUNDRY COLLECTORS, WITH OR WITHOUT INTERGRATED PUSH HANDLE

INTERGRATED PUSH HANDLE	
1 RING LAUNDRY COLLECTORS	6
2 RING LAUNDRY COLLECTORS	8
3 RING LAUNDRY COLLECTORS	10
4 RING LAUNDRY COLLECTORS	12
UNIVERSAL COMBINATION CART WITH DROPPY	17
OPEN LINEN WARD CARTS	25
LINEN WARD CARTS WITH CABINET	35
LINKABLE WARD CARTS	45
ACCESSORIES	53

ICONS

RVS Stainless steel frame

Swivel castors

Donut bumpers

Soft-close lid

Optional lid colour

Soot control

Push handle

S Linkable



1-RING LAUNDRY COLLECTORS

LC1W

1-RING LAUNDRY COLLECTOR / WALL MOUNTING

Stainless steel tubular frame, 22 mm (304).

Standard Care clamping ring.

On request also available with foot control.

ORDER NO.	WIDTH X DEPTH X HEIGHT
LC1W	390 X 340 X 290 MM





LC1

1-RING LAUNDRY COLLECTOR

Tubular stainless steel frame, 22 mm (304). Standard Care clamping ring.

ORDER NO.	WIDTH X DEPTH X HEIGHT
LC1	410 X 435 X 950 MM









1-RING LAUNDRY COLLECTOR

Tubular stainless steel frame, 22 mm (304).

Standard Care clamping ring.

Lid with damping & foot control.

Choice of lid colour (see page 55).

ORDER NO.	WIDTH X DEPTH X HEIGHT
LC1CDF	410 X 510 X 965 MM





1-RING LAUNDRY COLLECTOR

Tubular stainless steel frame, 22 mm (304). Standard Care clamping ring.

Includes donut bumpers & push handle.

ORDER NO.	WIDTH X DEPTH X HEIGHT
LCH1	490 X 455 X 950 MM



Detail of stainless steel tube with round weld.







LCH1CDF

1-RING LAUNDRY COLLECTOR

Tubular stainless steel frame, 22 mm (304).

Standard Care clamping ring. Includes donut bumpers & push handle, and lid with damping & foot control.

Choice of lid colour (see page 55).

ORDER NO.	WIDTH A DEPTH A HEIGHT
LCH1CDF	490 X 510 X 965 MM

2-RING LAUNDRY COLLECTORS

LC2

2-RING LAUNDRY COLLECTOR


Tubular stainless steel frame, 22 mm (304). Standard Care clamping ring.

ORDER NO.	WIDTH X DEPTH X HEIGHT
LC2	660 X 435 X 950 MM









3-RING LAUNDRY COLLECTORS















4-RING LAUNDRY COLLECTORS









The Care clamping ring securely mounts a wide variety of bag types.

Fold the ring down, hang the bag through the ring, and fold up the underside of the ring.

To remove the bag, simply fold the underside of the ring back down.







2-RING COMBINATION CART









3-RING COMBINATION CART









4-RING COMBINATION CART



















LOC22

LOC22

SNIJDERS CARE









4-RING WARD CART

mm. Platform spacing: 320 mm.

LOC24

Tubular stainless steel frame, 22 mm (304).

Two times two rings side by side, each with Care clamping

1410 X 600 X 1100 MM

ring. Three light grey platforms, each of size 560 x 480

OPEN LINEN WARD CARTS









OPEN LINEN WARD CARTS

















DISINFECTANT DISPENSERS

Hygienic dispensers for disinfectants, liquid soaps and lotions.

Capacities of 500 or 1000 ml. Dispense liquid manually or by automatic sensor technology (no contact).

Easy to clamp on all Snijders Care carts (stainless steel, 22 mm tube).

D500 Dispenser 500 ml

D500L Dispenser 500 ml with lock

D1000AST Dispenser 1000 ml with automatic sensor technology

*Optional: drip tray underneath dispenser.































LOC20LM1

LOC20LM1


SNIJDERS CARE

LINKABLE WARD CARTS


















LINKABLE LAUNDRY COLLECTORS









To connect: step on pedal - connect both carts - release pedal. To disconnect: step on pedal – disconnect carts – release pedal.

WD21L





LINKABLE DROPPY

For 575 x 1000 mm waste bags, with a stainless steel support $\,$ platform on the lower frame. Mounted on 75-mm swivel casters of which 2 Includes brake, with grey rubber profile, including donut bumpers at the corners.

ORDER NO.	WIDTH X DEPTH X HEIGHT
WD21L	550 X 590 X 980 MM









ACCESSORIES

LWR

CARE RING

This Care ring securely mounts a wide variety of bag types.

Fold the ring down, hang the bag through the ring, and fold up the underside of the ring.

To remove the bag, simply fold the underside of the ring back down.

LWR







DISINFECTANT DISPENSERS

Hygienic dispensers for disinfectants, liquid soaps and lotions.

Capacities of 500 or 1000 ml. Dispense liquid manually or by automatic sensor technology (no contact).

Easy to clamp on all Snijders Care carts (stainless steel, 22 mm tube)

D500 Dispenser 500ml

D500L Dispenser 500 ml with lock

D1000AST Dispenser 1000ml with automatic sensor technology

*Optional drip tray under dispenser.

D500 D500L







D1000AST



RECTANGULAR PLASTIC LIDS

Order no. of lid: CW1 - CG2 - CB3 - CR4 - CY5 - CG6 - CB7 - CB8

Order no. of lid with damping and hinge: CDW1 - CDG2 - CDB3 - CDR4 - CDY5 - CDG6 - CDB7 - CDB8



HINGE FOR COVER + D

Hinge for plastic lid.



FOOT PEDAL + F

Foot pedal.



Laurent Janssensstraat 105 5048 AR Tilburg, The Netherlands

T +31 (0)13 750 15 15

F +31 (0)13 463 86 35

info@snijderscare.com

www.snijderscare.com



